

# Thermomac

## Overhead Stirrer

### TM-II SERIES

Overhead stirrer is suitable for mixture of high viscosity liquid or solid-liquid. Widely used in chemical synthesis, pharmaceutical, physical and chemical analysis, petrochemical, cosmetics, health care, food, biotechnology and other fields.

#### Characteristics:



- High reliability, long continuous service, outstanding performance
- The speed is accurate and controlled, easy operation
- Brushless DC motor drive, low noise and maintenance free
- Overload protection, motor would automatically stop in case of continuous overload
- Anti-spill samples, PC control system starts up smoothly, prevent samples overflow
- CE marked
- 2 years warranty

## Technical Specification:

Model	TM-II-6C	TM-II-7C
Power(W)	60	80
Max. Capacity(L)	20	40
Voltage(V)	100-240	
Speed Accuracy(r.p.m)	±1	
Speed Display	LCD	
Speed Range(r.p.m)	100~2500	
Plate size(mm)	200*300	



## Ordering Information:

PART CODE	DESCRIPTION
001008003	TM-II-6C 20 Liters supplied with strap clamp, stand and shaft
001008004	TM-II-7C 40 Liters supplied with strap clamp, stand and shaft