

Thermomac

Refrigeration Capacity Recycable Coolers

TRC SERIES

Especially for small Rotary Evaporator with professional design and delicate figure which can effectively save space. It can be used on worktable or ground.

Characteristics:



- Circulating joint nozzle can rotate all round. Easy to connect with various system.
- All parts connecting with secondary refrigerant adopt stainless steel (SUS standard) and macromolecule anti-rot substance.
- Compressors of international famous brand with high reliability and can work for a long time.
- Combine with Rotary Evaporator and Vacuum Pump to form a complete small experimental
- 2 years warranty

Technical Specifications:

Model	TRC-400
Bath volume (L)	3
Temperature range (°C)	-15~25
Materials of stored liquid bath	304
Cooling capacity (W)	400
Flow (L/min)	17
Pressure (bar)	0.2
Overall Power (W)	500
Power Supply (V/Hz)	220/50
Dimensions (mm)	430L×245W×550H
Weight (kg)	26