

# Thermomac


## Water bath

## WB SERIES



This equipment applies to medical unit, college, scientific research institution and the department of chemical printing, dyeing, pharmacy, petroleum and electronic component production for the use of evaporation, arefaction, concentration and constant temperature heating.

## Characteristics:



- The inner chamber and covers are manufactured from stainless steel that processed against corrosion and leakage
- The outer surface is made of electrostatic powder steel.
- Intelligent temperature controller with PID program, digital tube display and function of timing and over-temperature protection.
- User friendly control panel including large bright LED displays for temperature and time
- Convenient drain for discharging the liquid in the tank quickly
- By the help of digital thermostat, a homogenous temperature distribution is obtained.
- Excellent temperature control of liquid for uniform and stable temperatures
- CE marked
- 2 years warranty

## Technical Specification:

MODEL	WB6	WB15	WB22
Capacity	6L	15L	22L
Tem. range	RT+5~100		
Tem. motion	$\pm 0.5^{\circ}\text{C}$		
Senility of Tem. Control	$\leq \pm 1^{\circ}\text{C}$		
Display error	$\leq \pm 2.5^{\circ}\text{C}$		
Inner Chamber size(w x d x h),mm	300 x 135 x 150	325 x 300 x 150	500 x 300 x 150
Exterior Size(w x d x h),mm	318 x 168 x 210	350 x 318 x 210	524 x 322 x 210
Net Weight(kg)	4.5	6	7.5
Power	500W	1000W	1500W
Power Supply	200-240V		



## Ordering Information:

PART CODE	DESCRIPTION
001015002	WB6 Water Bath 6.1 Liter Capacity
001015003	WB15 Water Bath 14.6 Liter Capacity
001015004	WB22 Water Bath 22.5 Liter Capacity