

Thermomac

DIGITAL ULTRASONIC CLEANER

TSC SERIES

Inner chamber, basket, top cover and housing are stainless steel.



DEGAS boosts ultrasonic power for faster removal of stubborn dirt.

Timer provides automatic shut-off at the set cleaning time.

SEMIWAVE lowers power to protect delicate items.

Can also be used as a water bath when ultrasonic mode is off.

Temperature and time can be set simultaneously from the large control panel.

Technical Specifications

Capacity	2 Litres	4.5 Litres	10 Litres	15 Litres	22 Litres	30 Litres
Ultrasonic Power	60W	180W	240W	360W	480W	600W
Heating Power	100W	200W	300W	400W	500W	600W
Thermostat	20 °C -80 °C	20 °C -80 °C	20 °C -80 °C	20 °C -80 °C	20 °C -80 °C	20 °C -80 °C
External Dimensions (mm)	175*160*210	325*180*225	325*265*280	360*325*285	530*325*285	530*325*325
Tank Dimensions (mm)	150*135*100	300*150*100	300*240*150	330*300*150	500*300*150	500*300*200
Operating Frequency	40 Khz	40 Khz	40 Khz	40 Khz	40 Khz	40 Khz
Inner Chamber	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Timer	0-30 min	0-30 min	0-30 min	0-30 min	0-30 min	0-30 min
Drain Valve	No	No	Yes	Yes	Yes	Yes
Gross Weight	2,9 Kg	5,5 Kg	8,4 Kg	10.1 Kg	13,5 Kg	15,41
Package Dimensions	260*260*290	410*255*310	435*360*350	460*410*340	635*405*340	635*405*390

