

Thermomac

Tablet Friability Tester

TAB10

Transparent drums are manufactured from high-quality acrylic. Samples are placed in the synchronised drums, which then operate at the selected settings.

MPU control enables precise adjustment of drum rotation speed and number of revolutions. The timer automatically stops the unit when the set time expires.



User-friendly, automatic and simplified operation provides accurate, precise and stable measurement performance.

Technical Specifications

Model	TAB10
Number of Drums	2 pcs
Drum Diameter	286 mm
Drum Depth	39 mm
Tablet Drop Height	156 mm
Rotational Speed	20-90 rpm
Speed Accuracy	±1 rpm
Number of Revolutions	10-900 rpm
Revolution Count Accuracy	±1 rpm
Power	220V/ 50Hz/ 60W
Dimensions	360 mm x 300 mm x 340 mm

