


Thermomac

Multi-Position Magnetic Stirring Heating Mantle

HM - HMS Series



Four or six positions enable simultaneous heating and stirring of multiple samples, saving time.



Technical Specifications

Model	HM450	HM401	HM650	HM601	HMS450	HMS401	HMS650	HMS601
Capacity	500 ml	1000 ml	500 ml	1000 ml	500 ml	1000 ml	500 ml	1000 ml
Operating Mode	Continuous							
Max. Surface Temperature	450 °C							
Temperature Control Range	RT+20~250 °C							
Heating Element	Nichrome Wire							
Stirring Range	-	-	-	-	0~1400 RPM	0~1400 RPM	0~1400 RPM	0~1400 RPM
Motor Power	-	-	-	-	0.01 kW	0.01 kW	0.01 kW	0.01 kW
Heating Power (kW)	4x0.25	4x0.35	6x0.25	6x0.35	4x0.25	4x0.35	6x0.25	6x0.35
Overall Dimensions (LxWxH)	690x180x180	810x210x200	1060x180x180	1240x210x200	690x180x180	810x210x200	1060x180x180	1240x210x200
Package Dimensions (LxWxH)	730x220x210	850x250x230	1100x220x210	1280x250x230	730x220x210	850x250x230	1100x220x210	1280x250x230
Voltage	200~240V 50/60Hz							
Net Weight (kg)	8	14	12	21	8	14	12	21
Gross Weight (kg)	10	16	15	24	10	16	15	24

