


Thermomac

Magnetic Stirrer Heating Mantle

HMS SERIES



Technical Specifications

Model	HMS250D	HMS500D	HMS001D	HMS002D
Capacity	250ml	500ml	1000ml	2000ml
Operating Mode	Continuous			
Max. Surface Temperature	450 °C			
Temperature Control Mode	PID Microprocessor			
Temperature Control Range	RT+20~250°C			
Speed Range	0~1600RPM			
Heating Element	Nichrome Wire			
Surface Material	High-Quality Aluminium Silicate			
Motor Power	0.01kW			
Heating Power (kW)	0.15	0.25	0.35	0.45
Overall Dimensions (W×D×H)	220x220x165	220x220x165	280x280x220	330x330x230
Package Dimensions (W×D×H)	230x230x195	230x230x195	290x290x300	340x340x350
Voltage	200~240V 50/60Hz			



Model	HMS250	HMS500	HMS001	HMS002
Capacity	250ml	500ml	1000ml	2000ml
Operating Mode	Continuous			
Max. Surface Temperature	450 °C			
Temperature Control Mode	-			
Temperature Control Range	-			
Speed Range	0~1600RPM			
Heating Element	Nichrome Wire			
Surface Material	High-Quality Aluminium Silicate			
Heating Power (kW)	0.15	0.25	0.35	0.45
Overall Dimensions (W×D×H)	200x200x160	200x200x160	260x260x200	260x260x200
Package Dimensions (W×D×H)	230x230x170	230x230x170	290x290x220	290x290x220
Voltage	200~240V 50/60Hz			

