

Thermomac

Hotplate Magnetic Stirrer TM Series



Technical Specifications

Model	TM14D
Maximum Temperature	320°C
Temperature Setting Range	5°C-320°C
Maximum Stirring Capacity	10 Litres
Speed Setting Range	200-2200 rpm
Speed Control	10 rpm / Incremental
Time Range	1-999 min / Continuous
Plate Material	Ceramic-Coated Stainless Steel
Plate Size	140 mm
Device Dimensions	248 x 156 x 115 mm
Power Output	600 W
Power Supply	110VAC, 60Hz / 230VAC, 50Hz

Accessories



4x20mm Blok



4x30mm Blok



7x16mm Blok



10x8mm Blok



12x4mm Blok



500ml R.B. Flask Block



Carrier Top Plate