

Thermomac

Hotplate Magnetic Stirrer

TM Series

Constant temperature, rapid heating and safety features, including automatic heater shut-off in case of flame.

Enables simultaneous heating and stirring safely and efficiently.

Stirring speed and heating temperature can be controlled and adjusted.

The TM12 features a simple control panel and is very easy to operate.

Heating power and stirring speed are adjustable in several steps. Maximum temperature: 380 °C; maximum speed: 2,000 rpm.



Technical Specifications

Model	TM12	TM45
Maximum Temperature	380 °C	320 °C
Speed Range	100-2000 RPM	-
Heating Plate Dimensions	120x120	350x450
Stirring Capacity	Up to 2,000 ml	-
Power	180W	1800W
Height with Support Rod	390 mm	-
Voltage	220/240V 50/60 Hz	220/240V 50/60 Hz
Dimensions (W x D x H)	200x120x90	400x500x220
Weight	1.9 Kg	30 Kg

