


Thermomac

Hotplate Magnetic Stirrer

TM Series



Technical Specifications

Model	TM18D
Maximum Temperature	550°C
Temperature Setting Range	5°C-550°C
Maximum Stirring Capacity	20 Litres
Speed Setting Range	200-2200 rpm
Speed Control	10 rpm / Incremental
Time Range	0-999 min / Continuous
Plate Material	Nano Crystal Ceramic
Plate Size	180x180 mm
Device Dimensions	313 x 205 x 110 mm
Power Output	1000 W
Power Supply	220-240V, 50Hz

