


Thermomac

Water Bath

WB Series



Technical Specifications

MODEL	WB6	WB15	WB22
Capacity	6.1 L	14.6 L	22.5 L
Temperature Range	RT +5 ~ 100 °C		
Temperature Variation	±0.5 °C		
Temperature Control Accuracy	≤ ±1 °C		
Display Error	≤ ±2.5 °C		
Inner Dimensions (W x D x H)	300x135x150	325x300x150	500x300x150
Outer Dimensions (W x D x H)	318x168x210	350x318x210	524x322x210
Net Weight (kg)	4.5	6	7.5
Power	500W	1000W	1500W
Power Supply	220/240V		

